



AFP Imaging Offers a Mobile CBCT Solution: The NewTom VG Mobile.

Some sensitive CBCT machines often require extensive calibration processes and timely warm-up periods making mobile applications difficult. With one button calibration and a simple, fast warm-up period, the NewTom VG makes going on the road easier.

NewTom VG Mobile can be configured to fit in a variety of van configurations to suit your needs. The NewTom VG Mobile from AFP Imaging is the CBCT machine that goes where you want your practice to go.

NewTom is the only system that will automatically adjust the radiation dosage based on patient size. This means that a small child should receive up to 40% less radiation than a full size adult. Other systems deliver the same amount of radiation regardless of patient size. With NewTom, radiation levels are different for every patient – low radiation and clear, accurate images.

Whether it is in your clinic or in your mobile office, AFP Imaging stands behind its products with a 2 year manufacturer's warranty, ensuring your practice goes wherever the road takes you.

For Your Patients
Significantly Lower Radiation

For Your Practice
The Highest Standard of Care

For You–NewTom VG Mobile
Full Service – One Unit

